

An *in vitro* T cell exhaustion model for characterization of multi-specific biologics

Lynne Dunsford¹, Rosie Thoires¹, Robert Benson¹, Thierry Giffon², Angus Sinclair², Agapitos Patakas¹.

¹ Antibody Analytics, Antibody Analytics Discovery Centre, Newhut Rd, Motherwell, ML1 3ST, UK. Email: apatakas@antibodyanalytics.com, ² IGM Biosciences Inc, Mountain View, CA

Introduction


Most studies investigating T cell exhaustion employ mouse models of chronic viral infections (e.g. LCMV), which are difficult or expensive to perform and, on most occasions, do not allow the assessment of human-targeting biologics.

Traditional T cell activation assays employing polyclonal stimuli (e.g. SEB-activation), do not generate exhausted T cells and are not suitable to assess whether a compound has an impact on T cell exhaustion.

There is an unmet need for a convenient biologically relevant *in vitro* model of human T cell exhaustion. Antibody Analytics have developed an *in vitro* model of chronic antigen stimulation that results in T cells with phenotypical and functional characteristics of exhaustion.


This model can be employed for the screening of immunomodulatory molecules, bispecific T cell engagers and multi-functional molecules.

Generation of Exhausted T cells




Isolated T cells were stimulated several times with antibodies targeting the TCR and costimulatory molecules. The generated phenotype of the exhausted T cells was assessed by flow cytometry and secondary functional assays

In vitro Tex cells express high levels of coinhibitory receptors




CD3+ T cells were either stimulated once or multiple times through their TCR and costimulatory molecules and the expression of several co-inhibitory and costimulatory molecules was assessed by flow cytometry. Similar results were observed in CD4+ and CD8+ isolated T cells and in multiple donors.

PD-1 blockade partially restores Tex cell functionality




Freshly isolated (non-exhausted) or exhausted T cells were used in a mixed lymphocyte reaction (MLR) in the presence of various concentrations of the PD-1 targeting antibody Nivolumab. IFN-gamma production was measured after 5 days by ELISA.

In vitro generated Tex cells are functionally impaired




In vitro generated Tex cells have diminished ability to produce IFN-gamma relative to non-exhausted (freshly isolated) T cells in a secondary T cell re-stimulation assay.



In vitro generated CD8+ Tex cells have diminished ability to induce target cell cytotoxicity mediated by a bispecific T cell engager.

A bifunctional antibody that targets PD-L1 and activates the IL-15 signaling pathway potently reverses T cell exhaustion



Treatment with IGM-7354 promoted increased secretion of IFN-gamma (A), proliferation (B) and T cell activation (C) by *in vitro* generated Tex cells in a MLR assay

Summary

- Robust generation of human *in vitro* Tex cells
- In vitro* generated Tex express high levels of co-inhibitory molecules
- In vitro* generated Tex cells are functionally impaired and have reduced ability to mediated effector functions
- PD-1 blockade partially restores functionality of Tex cells
- IGM-7354, an IgM-based bifunctional biologic, reverses T cell exhaustion more potently than anti-PD-L1 antibodies or a non-targeted IL-15 fusion molecule alone